



Low-Dead-Space Injection System

Proma inject Syringe + ZELOSTAT® 35G needle




Low-Dead-Space (LDS) Luer-Lock 1mL Proma inject Syringe AND 35G 4mm Low-Dead-Space Hypodermic Needle by ZELOSTAT® for Botulinum Toxin Micro Injections



MORE PRODUCT FOR YOUR PATIENT

The Low-Dead-Space technology based on a precise alignment between all components, dramatically reducing product waste.

From the first use you can already save up to 45µL botulinum toxin solution. Changing the needles 3 times will save over 100µL (0.1mL) of product equals to at least one extra injection point you can provide your patient.



KEEP YOUR ROUTINE

With an inner diameter of approx. 0.1mm, our 35G needle offers a flow rate and piston pressure similar to a conventional 34G needle.


Keep your injection technique — the injector will notice hardly any difference, the patient feels less needle.

SAFETY AND COMFORT

The Luer-Lock adapter tightly secures the injection needle on the syringe allowing the injector to work with a wider range of piston pressure without risking the needle to detach.

The specially designed plunger glides almost effortlessly, ensuring high level of control during the injection.

Our 35G needle is approximately 40% thinner* than conventional 34G needles and therefore can reduce the risk of injection-related side effects, thereby increasing patient comfort during the injection.



INNOVATION THAT PAYS

The efficiency in our state-of-the-art manufacturing facilities allows us to cut costs dramatically while maintaining the highest standards set by the American FDA and European EMA. We are proud to be able to forward the cost efficiency of our innovations to our customers offering the Proma inject Syringes at a price significantly below market average.

Visit us and order here:



The information provided in this data sheet is for general guidance only and may be subject to change. It is the user's responsibility to ensure the product is suitable for its intended purpose and used according to the manufacturer's instructions. Single use only. Do not reuse or re-sterilize.

PromaInject is part of PromaMedical GmbH
Wernher-von-Braun-Straße 10
85640 Putzbrunn – Munich,
Germany

T: +49 (0)89 7474713-0
F: +49 (0)89 7474713-1
M: info@promainject.com
W: www.promainject.com

e-mail us:



* 39,82% difference when comparing the cross-section area to TSK, THE INVISIBLE NEEDLE™